

# Junior Biologics formulation

Walk up, set up your run and walk away. Junior automates the manual measurements you do one by one at the bench. Choose a combo of high-throughput pH with visual inspection, visual inspection with viscosity or viscosity with pH. Crank through more samples in a single day, get them done the same way every time and check out a broader developability and formulation space.

### **Applications**

- Rapidly characterize a wide range of protein formulations with limited material
- Perform developability and preformulation screens
- Evaluate formulation robustness
- Formulation development
- Manage and track formulations and analytical results to facilitate rapid scientific decisions

### **Key features**

- Automate pH, viscosity and visual inspection
- Use a wide range of vials and microplates
- Check the pH of a 96-well plate in approximately 45 minutes
- Measure viscosity of protein formulations up to 100 cP with only 100 µL of sample
- Grab color, turbidity and visible particle count all in one shot
- Capture images and archive them for easy re-analysis



Junior configured for biologics formulation



Visual inspection station



Multi-channel pH probes



Viscosity station



Example Junior deck configured for biologics formulation

- 1 Vial/plate gripper
- 2 Waste bin
- 3 Viscosity station
- 4 Tool and plate rack
- 5 3-Position plate rack
- 6 Visual inspection station
- 7 Vial/plate hotel

# **Available options**

#### Visual inspection station (VIS) analyses

Includes:

· Visual particle analysis

Turbidity

Color measurement

Vial size: 2-20 mL

Recommended sample volume: 1 mL in 2 mL serum vial

Measurement time: 2-3 min per vial

#### Suspended visible particle detection

Minimum particle size detected: 80  $\mu m$  Maximum solution viscosity: 30–35 cP

Particle count accuracy:

No particles: 0 particles detected
1-3 particles: Detect at least 1 particle

• 4-9 particles: Actual particle count ±2 particles

• 10-25 particles: Actual particle count ±5 particles

#### **Turbidity**

Measurement range: 10-1000 NTU

Measurement accuracy:

O-100 NTU: ≤±5 NTU

>100-1,000 NTU: ≤5%

Repeatability: <1 NTU for 10 consecutive samples

#### Color measurement

Color: Correct match of Euro Pharmacopeia BY1-BY7 standards

#### pH measurement

Configuration: 4-channel glass probe

Measurement time per 96-well plate: <45 minutes

Range: 0-14 pH

Resolution: 0.01 pH unit Accuracy: ±0.03 pH unit Precision: ±0.05 pH unit

#### Viscosity station

Measurement range: 1-100 cP

Accuracy: ±0.5 cP + 10% of the actual viscosity Repeatability: StDev <0.5 cP + 5% of mean

Sample volume: 100 µL

Minimum volume in well: 200 µL Temperature range: 4–40 °C Temperature accuracy: ±1 °C Measurement time: 6 min/sample Throughput: 10 samples/h Vial/plate gripper

Plate size: Standard microtiter

Vial size: 1–125 mL Total mass: Up to 3 kg

Viscous liquid dispenser

Technology: Positive displacement tip (PDT)

Disposable tips:  $10 \mu L$  to  $10,000 \mu L$  from Eppendorf and Rainin

Viscosity: 1 cP to 1,000 cP

Vortexing station

Orbital: 60-3570 rpm

Maximum vortexing mass: 2268 g (5 lb/plate)

#### Off-deck third-party instrument virtual integration

• DLS

HPLC

cIEF

Other systems available for virtual integration.

Please contact Unchained Labs for a full list of systems.

# **Facilities requirements**

#### Physical

• With integrated enclosure: 105 cm W x 90.4 cm D x 140 cm H ~150 kg

With integrated table option:

167 cm W x 90.4 cm D x 200 cm H ~240 kg

#### **Electrical**

• Junior:

110-220 ±10% VAC, 50-60 Hz, 16 A

Computer:

US: 115 V ±10 %, 60 Hz, 10 A

EU: 220-230 V ±10 %, 50 Hz, 16 A Compressed dry air: 0.55-0.9 MPa (80-130 psi), 40 L/s

**L**ABS

6870 Koll Center Parkway

**Unchained Labs** 

Pleasanton, CA 94566 Phone: 1.925.587.9800 Toll-free: 1.800.815.6384 Email: info@unchainedlabs.com